




Toshiro Shinke<sup>3</sup>, Liam P. Mullins<sup>1</sup>, Norah O'Brien<sup>1</sup>, Kevin B. Heraty<sup>1</sup>, Tony McMahon<sup>1</sup>, Martin G. Burke<sup>1</sup>, Chas Taylor<sup>1</sup>, Paul Gilson<sup>1</sup>, Keith Robinson<sup>3</sup>, Colin G. Caro<sup>1,2</sup>,  
1, Veryan Medical Ltd, Units 5/17 City Business Centre, Brighton Rd., Horsham, West Sussex, RH13 5BB, UK  
2, Dept. of Bioengineering, Imperial College London, London SW7 2AZ, UK 3, Saint Joseph's Translational Research Institute, Atlanta, GA, USA

## Background

- Swirling blood flow patterns in the vasculature augments vessel wall shear stress, which in turn may have a profound effect to suppress neointima formation post stent implant.
- Veryan Medical, Ltd has developed a novel peripheral stent with 3D centreline curvature

## Objective

- The objectives of this study were to determine, using a porcine model of carotid stenting:
  - Evaluation of BioMimics 3D stent deliverability to assess any issues at time of implant
  - Evaluation of angiographic and duplex ultrasound data to assess a trend towards an increase in swirling flow in BioMimics 3D stents relative to controls at implant and pretermination
  - Evaluation of histological endpoints to assess if there is a trend toward lower neointimal thickness and % stenosis in the BioMimics 3D stents




Example of a stent with BioMimics 3D technology

## Methods

- Test system: Juvenile domestic Yorkshire hybrid pigs (N=10, body weight 42.7±3.5 Kg)
- Stent implant procedure:
  - Animals were implanted with one BioMimics 3D stent in a carotid artery and one straight stent in contra-lateral carotid artery according to standard technique
  - Mean lumen diameter (mLD) within the stent were measured
- Angiographic restudy, termination, & specimen harvest:
  - Angiographic restudy at one-month post-implant
  - After termination, the carotid arteries were excised and processed for histological analysis
- Color Doppler and Duplex ultrasound examination
  - Assess a trend towards an increase in swirling flow in the BioMimics 3D stents relative to control segment
- Stented carotid segments were embedded in methyl methacrylate; histologic sections were prepared using the precision saw-and-grinding method
- Histologic sections were evaluated by computer-assisted planimetric histomorphometry to measure:
  - area of arterial lumen
  - internal elastic lamina (IEL) area
  - thickness of neointima
  - area of neointima (IEL area – lumen area)
  - % area stenosis  $[(1 - (\text{luminal area} / \text{IEL area})) \times 100]$
- Histomorphologic assessment of vessel wall response


## Results

### ➤ Angiography Immediately Post-Implant




Swirling blood flow pattern was observed in a BioMimics 3D stent immediately after implantation, using small-volume contrast injection. In straight stents, swirling flow pattern was not detected.

### ➤ Duplex Imaging Post-Implant




Swirling blood flow pattern was evident by the presence of bi-directional flow in a BioMimics 3D stent immediately after implantation. Orientation of flow is uni-directional i.e. no swirl is present in straight stents.

### ➤ Angiography One Month Post-Implant




mLD was similar between the groups after implant (HS: 5.0±0.4mm, SS: 5.0±0.4mm; P=0.903), while all animals showed higher mLD in HS at one-month (HS: 5.2±0.6mm, SS: 4.7±0.6mm; P=0.053).


### ➤ Histology of Carotid Stents at One Month Post-Implant



### ➤ Histomorphometry



### ➤ Histopathologic scoring



- Neointima within stented segments consisted of smooth muscles in an organised extra-cellular matrix, with uniform endothelial cell coverage both in BioMimics and Straight stent.
- Inflammation score was equivalent between BioMimics 3D and Straight stent

## Summary

- BioMimics 3D stents were readily deployed in straightforward manner using standard equipment and supplies, into pig carotid arteries
- BioMimics 3D stents were capable of imposing a helical profile on a porcine disease free common carotid artery
- The subsequent reshaping of the vessel induces swirling flow as verified using duplex Doppler ultrasound and digital subtraction angiography
- Histopathology studies indicate that the helical stents reduces levels of neo-intimal thickness by 45% (P<0.001) and stenosis by 34% (P<0.001) when compared to straight stents in a porcine carotid model.
- Extensive endothelialisation of the stented vessel segment has taken place at one-month

## Conclusions

- BioMimic 3D stent provides a novel potential therapeutic strategy for vascular intervention
- BioMimics 3D design promotes swirling blood flow and inhibits neointima formation in stented porcine carotid arteries, possibly due to the enhancement of vessel wall shear stress
- Neointimal morphology was mostly similar to standard stent